



FLEXIBILITY IN A WORLD OF MIXED ABILITY

Overnight Instructions



DO NOT RE-CONFIGURE THE CHAIR BEFORE BEING SHOWN HOW TO DO SO.



KEEP HANDS AWAY FROM THE EDGES OF THE CUSHIONS WHEN RE-CONFIGURING THE CHAIR.



Bed to Chair



Lift up front edge of seat with handle.




Pull seat unit forward letting it roll on floor. No need to lift!



When fully extended the bed surface is flat and level.

Chair to Bed



Push down in in the middle of the edge of the main cushion.



Push the middle section to form the 'back' of the chair



Push the chair back until the front of the seat drops into position.

Sleep with feet towards the handle on the front of the chair.

Do not sit on the 'head' end of the bed.

Do not lift the seat section off the floor at any time!

For cleaning and decontamination procedures, refer to the appropriate section of the Fabric and Cleaning Guide supplied with the chair. Careflex use a range of covering materials, each with particular characteristics and specific cleaning guidelines. If you are unsure of the chair's covering material, contact CareFlex with the serial number for the chair which will allow CareFlex to trace the specification of the chair.