CareFlex

Product Storage Guide UI

PRODUCT STORAGE SUPPLEMENT

WARNING

If you are in any doubt as to the safe use of your CareFlex product, or if any parts have failed, please cease using the product and contact CareFlex or your local dealer as soon as possible



AlphaCare II



EasiTilt



HydroCare



HydroTilt



HydroTilt XL



MultiAdjust



MultiAdjust Advance



HydroFlex



SmartSeatPro II



Important Information

- CareFlex products are designed for indoor use and when not in use should be stored in a dry place that is not subject to extremes of temperature.
- Clean the product in line with the instructions in the Fabric & Cleaning Guide before storage/re-issue.
- Carry out basic maintenance checks to ensure your product is in good working order prior to storage or re-issue, checking function and condition of:
    - Castors and brakes.
    - Adjustment of TiS/BAR/Seat-depth/Arm width and back/wing articulation, as appropriate to specific product.
    - Cushions/covers, checking for rips & tears.
    - Check chair runs straight when pushed.
    - Visually check the structure of the product, paying particular attention to weld points on the frame, ensuring there are no signs of fatigue or cracking around the welds.
- When the product is stored or left stationary, ensure that all castors are locked (if braked) and facing away from the base/chassis, as this improves product stability. This is especially important when Tilt-in-Space or Back Angle Recline adjustment has been made.
- Never leave a product on a sloping or uneven surface, or in an area in which it would provide obstruction to others.
- Powered product can be charged using the supplied adaptor. Do not leave the product plugged in once fully charged, and do not position the product as to cause an obstruction while charging.
- CareFlex products comply with fire safety regulations. However, it is best practice to ensure the product is kept away from all direct sources of heat including naked flames, cigarettes and gas heaters.



If you are in any doubt as to the safe use of your CareFlex product, or if any parts have failed, please cease using the product and contact CareFlex or your local dealer as soon as possible

CareFlex Ltd.
Templer House, King Charles Business Park,
Old Newton Road, Heathfield,
Devon TQ12 6UT

Tel: 0800 018 6440 | www.careflex.co.uk | info@careflex.co.uk